




Plant Gene Editing with Cas-CLOVER™
“the clean alternative to CRISPR/Cas9”

What is Cas-CLOVER?

Hybrid gene editing system using nuclease-inactivated Cas9 protein fused to Clo51 endonuclease.

dCas9 serves as a fusion protein only. It is mutated and unable to cut DNA.

Cleavage activity depends on the dimerization of an “obligate dimer” of Clo51.



Benefits



Ease of use: double gRNA guided

High efficiency: dCas9-guided RNA recognition


High fidelity: cuts only when Clo51 nuclease dimerizes

Optimal Activity Reported in ~16-30 Base Pair Range


Designing multiple Cas-CLOVER sites within target genes is straightforward.



Extensive studies reveal minimal off-target activity



- Left y-axis = percent knockout (indels) for the given gRNA combination
- Right y-axis = number of next-generation sequencing (NGS) reads. Comprehensive coverage from around 100,000 reads
- X-axis bar graphs = percent knockout (indels) the different gRNA combinations for predicted off-target and on target "ON"
- X-axis lines = dotted blue line is 1 standard deviation (SD) from background mutation rate. Solid blue line is 2 SD from background mutation rate



- Other studies examine few (10 to 25) candidate off-target sites¹⁻³.
- Our Cas-CLOVER off-target study is ~10x broader and includes 8 donor lots.

Cas-CLOVER plant genome engineering



Robust genetic
tools for trait
development

Clean editing for
possible non-GMO
crop designation

Optimization of the gene editing system:
targeting the PDS gene in tobacco for visual
phenotype read-out

Utilized visual PDS read-out to optimize cutting
efficiencies to 100%, demonstrating editing of
all 4 alleles

T0 and T1 generations were screened for
stable gene editing

First proof-of-concept: PDS Knockout


Knockout mutation of the phytoene desaturase (PDS) gene results in a visually pale and white phenotype



Green, pale mixed, and white shoots were all screened to determine cutting efficiencies. Afterwards, we verified that only white and pale material was edited at an overall average of around **50%+ editing efficiency (white shoots/total shoots from plates)**.

Optimizing cutting efficiencies of Cas-CLOVER in plants


Demeetra developed an optimized protocol which is transferred to our licensees and partners



Improved system consistently yields 100% editing efficiency, modifying all four alleles in tobacco

Editing Over Time

Percent of total population edited over time




- Correctly targeted plants can be identified in a relatively short period of time
- Editing efficiencies plateaued at the Trim 3 timepoint

Green = no editing
Yellow = 1-2 alleles edited
White = 3-4 alleles edited
Note: Yellow plants can be crossed to produce fully edited plants.
Therefore, the data represent white phenotype as “% fully edited” and white + yellow “% edited”

Whole callus editing identification

Editing % of PDS in whole callus (C1-7)



Next, to speed up and simplify our process, we aimed to test whole callus from plates for cutting efficiency with Cas-CLOVER

- Whole callus were homogenized, DNA was extracted, PCR amplified, sequenced, then analyzed by editing software

- These data are from a single promoter expressing Cas-CLOVER, Demeetra has tested multiple promoters with varying efficiencies.



Examples of edited sequences

POWERED BY SYNTHEGO ICE


RELATIVE CONTRIBUTION OF EACH SEQUENCE (NORMALIZED)



Relatively large deletions (in comparison, CRISPR/Cas9 generates 1-3 bp deletions) is consistently observed with Cas-CLOVER

This results in simpler, faster and more cost-effective screening of plants for editing, see [our post](#) on this.

Stability of plant edits over time




- Red annotated plants were propagated to maturity
- Mature plants in soil demonstrated stable editing levels at flowering stage

T1 generation retain Cas-CLOVER edits

- Selected plants were self pollinated, and seed was harvested
- Unpurified PCR products screened on QIAxcel for rapid PCR detection

QIAxcel T1 Plant Knockout Gel



QIAxcel Mutagenesis Report

Individual Plant Samples	Deletion reads						WT reads			
	Pos	Sample Info	Found	Size [bp]	Conc. [ng/μl]	Rel. Abund. [%]	Found	Size [bp]	Conc. [ng/μl]	Rel. Abund. [%]
	E2	21B2.2	Yes	440	0.49	13.55	No	n/a	n/a	n/a
	E4	21B2.5	Yes	438	0.97	26.83	Yes	468	0.37	7.46
	A2	21B2.9	Yes	448	0.18	4.83	Yes	456	0.17	3.45
	A3	21B2.10	Yes	437	0.46	12.73	No	n/a	n/a	n/a
	A4	21B2.11	Yes	446	2.66	73.35	No	n/a	n/a	n/a
	A5	21B2.12	Yes	439	1.65	45.53	No	n/a	n/a	n/a
	A6	21B2.13	Yes	441	0.56	15.59	No	n/a	n/a	n/a
	D11	21B2.14	Yes	443	0.18	4.89	Yes	463	0.39	7.80
	A7	21B2.15	Yes	442	1.31	36.18	No	n/a	n/a	n/a
	A8	21B2.16	Yes	448	1.51	41.64	No	n/a	n/a	n/a
	B5	21B3.10	No	n/a	n/a	n/a	No	n/a	n/a	n/a
	B9	21B4.1	Yes	438	0.54	14.83	No	n/a	n/a	n/a
	D8	21B4.8	No	n/a	n/a	n/a	Yes	469	0.08	1.67
	C10	21B8.10	Yes	445	2.11	58.27	No	n/a	n/a	n/a
	C11	21B8.11	Yes	444	1.80	49.58	No	n/a	n/a	n/a
	D1	21B8.13	No	n/a	n/a	n/a	Yes	473	4.98	100.0
	D3	21B8.15	Yes	444	2.05	56.63	No	n/a	n/a	n/a

- Cas-CLOVER produces larger indels than other technologies
- The large indels are rapidly and easily detected by peak size & translated into a gel image

Deletions double checked by Synthego ICE

T1 generation screened by ICE

- 14 of 16 samples sequenced and analyzed demonstrated 100% indel percentage
- the two noted samples below also showed detection of WT remaining by QIAxcel
- two other samples with noted WT peaks by QIAxcel revealed only very small deletions by sequence



Licensing Options

Demeetra AgBio holds the exclusive license in Agriculture to an issued patent portfolio for Cas-CLOVER. We offer sublicenses and strategic development partnerships.



Evaluation

Technology transfer with reagents and protocols for a 6-month evaluation research license

Research

Economically attractive yearly fee. Freedom to switch technologies for commercialization to navigate gene editing regulatory landscape

Commercial

Flexibility of one-time fee-based license or more standard upfront/milestone/royalty structure

For More Information, Contact



Jack Crawford, M.S.

CEO, Demeetra AgBio, Inc.

Mobile: 315-351-9115

Email: jcrawford@demeetra.com

