



Demeetra


**Complementary Gene Editing Technologies Enhance
Rapid Generation of Stable Cell Lines Producing
Enveloped Virus-Like Particles**

Corey Brizzee, PhD
Director of Gene Editing

Demeetra's Core Technologies




piggyBac Transposase/transposon



- Transposase integrates cargo at highly expressed sites – open chromatin and non-coding regions
- Used as a genome surveyor tool

Cas-CLOVER Targeted Nuclease



- High efficiency targeted dimeric nuclease
- High fidelity: cuts only when Clo051 nuclease dimerizes
- Larger deletions with overhanging staggered cuts – high knock-in efficiency


Bioproduction Platforms for Rapid Response to Pandemics

Overall Goals:

- 1) Stable Cell Line expressing Core GAG protein
 - Increased levels of GAG correlate to higher eVLP production
- 2) Rapid **and** Tailored expression of Antigenic Proteins
 - Specific ratios of expression lead to better efficacy in animal studies


To achieve these goals:

- 1) piggyBac transposase to identify novel sites in genome for a wide range of GAG expression
- 2) Targeted knock in of VLP cargo elucidates site-specific expression contribution of novel sites compared to known genomic safe harbors





A perfect match for enveloped viruses

Novel High-Expressing Site Identified via piggyBac Genome Mapping



HEK293SF Single-Cell Clonal Population



TOPO Clone
↓
Colony PCR/Sequence

Clone Site Mapping and GAG expression

Clone #	Total Integrations	[GAG] (ng/ml)	Genomic Locations
1 (b2)	4	0.0629	Chr.8 (73,738,168)
3 (b4)	4	0.0670	Chr.10 (11,461,767)
4 (b5)	4	0.1596	Chr.17 (21,100,584)
6 (b7)	4	0.1238	Chr.22 (41,468,704)
7 (b8)	4	0.0373	
8 (b9)	1	0.1001	Chr.2 (proprietary)

Cas-CLOVER Targeted Knock-In Concept and Workflow


Novel site in Chromosome 2 was compared to known Safe-Harbor sites, **Rogi1**(high) and **GSH31**(medium), from Aznauryan et al. 2022

- Typically design three guide pairs per target region
- Test with various ratios of guides to Cas-CLOVER

Cell Reports Methods


Discovery and validation of human genomic safe harbor sites for gene and cell therapies

Erik Aznauryan,^{1,2,3,4} Alexander Yermanos,^{1,5,6} Elvira Kinzina,⁷ Anna Devaux,⁸ Edo Kapetanovic,¹ Denitsa Milanova,^{3,4} George M. Church,^{3,4,9} and Sai T. Reddy^{1,9,10,*}




At Demeetra, we simultaneously test guide while performing all knock-ins to save 5-7 days.

Site-Specific Knock In of eVLP Cargo Validates Novel Chromosome 2 Site



Adapted from Aznauryan et al. 2022, *Cell Reports Methods*



actatagggcga
tgatatcccgct
ter

Sample Identification	VCD (E+6 cells/mL)	Viability (%)	[Gag] (ng/mL)
GSH31	2.08	95	3724
Chromosome 2	1.58	97	3535
Rogi1	1.83	92	<LOD at 1/100

actatagggcgaattggagctccaccgCGGTGGCGCCGCTCTAGAGAGCTTGGCCCATTCATACGTTGTATCCATATCATAAATATGTACAT
ACTATAGGGCGAATTGGAGCTCCACC GCGGTGGCGCCGCTCTAGAGAGCTTGGCCCATTCATACGTTGTATCCATATCATAAATATGTACAT

Original Sequence
2-4_0.3P_NF072523_5Ta

ACCTTATAAATAATCATTCTGTTTTATTCTTAAATCCCTACAGTACACTGAGCATATAAACTAAGGAAATTTCCATGAGCATAATAAGGAGATTATTTATCTGCTTAAATCTTTCAGTTATTCTGTGATGCTC'
TGGAATATTTATTAGTAAGACAAAAATAAGAATTTAGGGATGTCATGTGACTCGTATATTTGATTCCTTTAAAGGTA CTGATTTATCTCTAATAAATAGACGAATTTAGAAAGTCAATAAGACACTACGAG,
3' ARM OF HOMOLOGY
Original Sequence
2-4_0.3P_NF072523_3Targ_nRogi3trgste2-R4_G54433_4
ACCTTATAAATAATCATTCTGTTTTATTCTTAAATCCCTACAGTACACTGAGCATATAAACTAAGGAAATTTCCATGAGCATAATAAGGAGATTATTTATCTGCTTAAATCTTTCAGTTATTCTGTGATGCTC'

Current Work and Future Platforms

Develop Stable Cell Lines in a timely manner


- Leverage strategies to allow for rapid integration of newly identified future variances of concern, or “Disease X”.
- Newly developed Bivalent cell lines produced in 34 days
- Future bicistronic and polycistronic knock-ins for new product candidates for manufacturing novel vaccine therapies.

Multiple sites or single site
Tailored ratios


Quadrivalent




Stable GAG Cell Line




Trivalent



Bivalent



Monovalent



Questions?