

CleanCut GS CHO Single-cell Cloning SOP



SOP003 | Version: 1E01 | Effective Date: 08/20/2026

Purpose

This SOP outlines a general workflow for the preparation and handling of CleanCut GS CHO cells before and after single-cell cloning. It is intended to support a variety of single-cell cloning methods and platforms, including limiting dilution cloning and cell-dispenser instrument-based workflows. Detailed procedures for performing specific single-cell cloning methods are outside the scope of this SOP.

Required Materials

Product or Reagent	Vendor	Catalog No.
Fisherbrand™ Disposable PES Filter Units	Fisher Sci	FB12566506
GlutaMAX™ Supplement	Gibco	35050061
CD FortiCHO™ Medium	Gibco	A1148301
Anti-Clumping Agent	Gibco	0010057AE
HyClone™ Penicillin-Streptomycin 100X solution	Cytiva	SV30010
125 mL Polycarbonate Erlenmeyer Flask (E-125)	Corning	431143
Solentim® InstiGRO CHO PLUS	Advanced Instruments	RS-1225
Nunc™ Non-Treated 24-Well Plate	Thermo Sci	12-566-82
25 cm ² Rectangular Canted Neck Not Treated Cell Culture Flask (T-25)	Corning	431463

Media Preparation

Complete CHO Medium (1L)

Add the following components to a 0.2 µm aPES 1L filter unit. Filter by suction once all components are added.

- 40mL GlutaMAX™ Supplement (final conc – 8mM)
- 5mL Penicillin-Streptomycin (*Note: this component is optional*; final conc – 0.5%)
- 1mL Anti-Clumping Agent (final conc – 1:1000)
- Bring total volume to 1000mL with CD FortiCHO Medium

GS Selection Medium (1L)

Add the following components to a 0.2 µm aPES 1L filter unit. Filter by suction once all components are added.

- 1mL Anti-Clumping Agent (final conc – 1:1000)
- 5mL Penicillin-Streptomycin (*Note: this component is optional; final conc – 0.5%*)
- Bring total volume to 1000mL with CD FortiCHO Medium

❖ *Note: For cells that underwent GS selection, we recommend continuing to use GS Selection Medium for single-cell cloning.*

Single-Cell Cloning Medium (50mL)

Combine InstiGRO CHO PLUS with either Complete CHO Medium or GS Selection Medium to a final concentration of 5%.

- Add 2.5mL InstiGRO CHO to the bottom of a 50mL conical (final conc – 5%)
- Add 47.5mL media.
- Vortex to mix.

❖ *Note: From this part of the protocol forward, any mention of “media” that is not Single-Cell Cloning Medium refers to either Complete CHO Medium or GS Selection Medium.*

Standard CleanCut GS CHO Cell Incubation Conditions

Temperature: 37°C

CO₂ level: 8%

Shaking:

- Well plates – Do not shake
- T-25 Flasks – 90 RPM, 19mm Orbit
- E-125 Flasks – 125 RPM, 19mm Orbit

Protocol

1. One day prior to single-cell seeding, passage cells to 1x10⁶ cells/mL in 30mL in an E-125.
2. On the day of seeding, use the single-cell seeding method of choice to seed 1 cell in 200µL of Single-Cell Cloning Medium per well into a 96-well plate.

3. For monoclonality confirmation, image the wells within 5hr of seeding.
 4. Allow cells to incubate for 14 days, imaging the wells again on day 7 to monitor growth.
 5. On Day 14, image the cells again and feed with 60 μ L of fresh media.
 6. On Day 15, consolidate clones with verified monoclonality and confluency of at least 33% to a new 96-well plate. To do this, gently mix each well and transfer the entire volume to the new plate.
 7. On day 17, transfer each clone from the 96-well plate to a low-attachment 24-well plate. To do so, add 400 μ L media to each well of the new plate and transfer the full volume of cells.
- ❖ *Note: Alternatively, add 200 μ L of fresh media to the 96-well plate, mix, and transfer 200 μ L of cells to a new 96-well plate. Use the cells from one plate for DNA extraction for genotyping and transfer corresponding clones with correct genotype to a 24-well plate 24-72hr later.*
8. When total cell count reaches at least 1×10^6 cells, transfer cells into a T-25 with 5mL media. Shake at 90 RPM.
 9. When Cell number reaches at least 3×10^6 cells, transfer into E-125 shaker flask 0.3×10^6 /mL in 10-30 mL media on shaker at 125 RPM.



Demeetra AgBio, Inc.

2277 Thunderstick Dr., #300, Lexington, KY 40505

For more Demeetra resources:

<https://demeetra.com/resources/guides-protocols/>

Contact us:

<https://demeetra.com/contact-us/>
